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Aim

Comparison of different antigen rapid tests using uniform sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs which are used as test specimens included in an evaluation panel.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable Cq values diluted 1:10 in negative samples in PBS (Cq = Quantification cycle; also reported as CT = cycle threshold). The Cq values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standard. In the case of the PCRs used, a Cq value of 25 corresponds to around 10^6 RNA copies/mL. The replication of the virus in cell culture as another characteristic of the samples was determined as a possible correlate for infectiousness.

Due to the high number of tests under evaluation few specimens of the first panel were replaced (resulting in panel 1 versions 1 and 2), and a new panel (panel 2) was manufactured in analogous manner. For panel 1 versions 1 and 2, 18 or 17 samples, respectively, were characterized as Cq ≤ 25 , 23 samples as Cq 25-30, and 9 or 10 samples, respectively, as Cq ≥ 30 ; the composition of panel 2 is: 20 samples with Cq ≤ 25 , 20 samples with Cq 25-30, and 10 samples with Cq ≥ 30 .

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 μ l of the pool (panel 1 versions 1 and 2), or 25 μ l (panel 2) were analysed using the components of the test provided, e.g. swabs and buffer. In case of tests designed for various sample types (e.g. nasopharyngeal, nasal, saliva, etc), the test specific swab dedicated for nasopharyngeal samples was used.

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Kontakt:

E-Mail: sarscov2ivd@pei.de

Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers using the uniform sample sets allows an overview of the current state of the art regarding sensitivity.

The current state of the art was defined as corresponding to a minimal sensitivity of 75% for the pools with $Cq \leq 25$. The results do not allow any conclusions regarding specificity of the tests.

Table 1 summarizes the results of those POCTs assessed as reflecting the current state of the art (alphabetical order of manufacturers).

Table 2 summarizes the results of those POCT (alphabetical order of manufacturers) which do not fulfill the sensitivity criteria and consequently are deleted from the list of reimbursable tests (see the homepage of Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM).

The first part of this evaluation has already been published (1); this publication summarizes further details of the comparative evaluation.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for reimbursement, and that few other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

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